












## 12.1. Bondable Buccal Tubes



### UPPER Tooth

	Prescription	Torque	Angle	Distal Offset	Length (mm)	Occlusal or Gingival	Hdgr Tube	Reference Code			
								Slot 0,18		Slot 0,22	
								Right	Left	Right	Left
1° Molar Upper-Triple Convertible with Hook NO Hook	 Roth	-10	0	+8°	3.6	O	.045	351-830	351-831	351-936	351-937
	MBT	-14	0	+10°	3.6	O	.045	351-865	351-866	351-965	351-966
	Ricketts	0	0	0	3.6	G	.045	351-011	351-012		
	Ricketts	0	0	0	3.6	G	.045	351-013	351-013		
1° Molar Upper - Double Convertible whit Hook	 Roth	-10	0	+8°	3.6	G	.045	351-823	351-824	351-901	351-902
	MBT	-10	0	+8°	3.6	O	.045	351-832	351-833	351-932	351-933
	Edge STD	-14	0	+10°	3.6	O	.045	351-834	351-835	351-934	351-935
Ricketts	 Ricketts	0	0	0	3.6	O	.045	351-805	351-806	351-905	351-906
1° Molar Upper-Single Convertible Whit Hook	 Ricketts	0	0	0	3.6	O	.045	351-015	351-015	351-115	351-115
	Roth	-10	0	+8°	3.6	-	-	351-821	351-822	351-921	351-922
	MBT	-14	0	+10°	3.6	-	-	351-873	351-874	351-973	351-974
	Edge STD	0	0	0	3.6	-	-	351-007	351-008	351-107	351-108
1° Molar Upper-Single Not-Convertible Whit Hook	Bi-Dim	0	0	+8°	3.6			351-005	351-006	351-105	351-106
	Roth	-10	0	+8°	3.6	-	-	351-215	351-216	351-315	351-316
	MBT	-14	0	+10°	3.6	-	-	351-273	351-274	351-381	351-382
	Edge STD	0	0	0	3.6	-	-	351-004	351-005	351-104	351-105
Twin Bracket conv.	Bi-Dim	0	0	+8°	3.6			351-001	351-002	351-101	351-102
	 Edge STD	0	0	0	3.6	-	-	301-201	301-201	301-301	301-301

### LOWER Tooth

	Prescription	Torque	Angle	Distal Offset	Length (mm)	Occlusal or Gingival	Hdgr Tube	Reference Code			
								Slot 0,18		Slot 0,22	
								Right	Left	Right	Left
1° Molar Lower Double Rect.	 Roth-Rect.	-25	0	5°	3.8	G	-	351-828	351-829	351-928	351-929
	MBT-Rect.	-20	0	0	3.8	G	-	351-819	351-820	351-923	351-924
1° Molar Lower-Double Round Conv. whit Hook	 Roth-Round	-25	0	5°	3.8	G	.045	351-817	351-818	351-917	351-918
	MBT-Round	-20	0	0	3.8	G	.045	351-291	351-292	351-391	351-392
Ricketts	 Ricketts	0	0	0	3.8	O	.045	351-015	351-015	351-115	351-115
1° Molar Lower Single Convertible Whit Hook	 Roth	-25	0	5	3.8	-	-	351-813	351-814	351-913	351-914
	MBT	-20	0	0	3.8	-	-	351-815	351-816	351-915	351-916
	Bi-Dim	0	5°	0	3.8	-	-	351-007	351-008	351-107	351-108
	Edg Std.	0	0	0	3.8	-	-	351-005	351-006	351-105	351-106
2° Molar Lower-Single Not Convertible with Hook	 Roth	-25	0	5°	3.8	-	-	351-653	351-654	351-335	351-336
	MBT	-10	0	0	3.8	-	-	351-293	351-294	351-393	351-394
	Bi-Dim	0	0	5°	3.8			351-003	351-004	351-103	351-104
	Edg Std.	0	0	0	3.8	-	-	351-001	351-002	351-101	351-102

## 12.2. Bondable Tubes Anatomic Base

### 1°/2° Molar Single Not Convertible w/h



Maxillary	Torque	Angle	Offset	Tube	Reference Code		Reference Code	
					Slot 0,18		Slot 0,22	
					Right	Left	Right	Left
Edgewise Standard	0	0	0	3.6	352-001	352-002	352-101	352-102
Bi-Dimensional	0	0	5°	3.6	352-003	352-004	352-103	352-104
MBT	-14°	0	10°	3.6	352-281	352-282	352-381	352-382
Roth	-10°	0	0	3.6	352-277	352-278	352-377	352-378
Roth	-10°	0	8°	3.6	352-215	352-216	352-315	352-316

Mandibular	Torque	Angle	Offset	Tube	Reference Code		Reference Code	
					Slot 0,18		Slot 0,22	
					Right	Left	Right	Left
Edgewise Standard	0	0	0	3.6	352-007	352-008	352-107	352-108
Bi-Dimensional	0	0	5°	3.6	352-009	352-010	352-109	352-110
Bi-Dimensional	0	0	8°	3.6	352-011	352-012	352-111	352-112
MBT 1° Bicuspid	-20°	0	0	3.6	352-295	352-296	352-395	352-396
MBT 2° Bicuspid	-10°	0	0	3.6	352-293	352-294	352-393	352-394
Roth	-25°	5	0	3.6	352-233	352-234	352-335	352-336

## 12.3. Bondable Tubes Big Fit Base

### 1°/2° Molar Single Not Convertible w/h








Maxillary	Torque	Angle	Offset	Tube	Reference Code		Reference Code	
					Slot 0,18		Slot 0,22	
					Right	Left	Right	Left
Edgewise Standard	0	0	0	3.6	353-001	353-002	353-101	353-102
Bi-Dimensional	0	0	5°	3.6	353-003	353-004	353-103	353-104
Bi-Dimensional	0	0	8°	3.6	353-005	353-006	353-105	353-106
MBT	-14°	0	10°	3.6	353-281	353-282	353-381	353-382
Roth	-10°	0	0	3.6	353-277	353-278	353-377	353-378
Roth	-10°	0	8°	3.6	353-215	353-216	353-315	353-316






Mandibular	Torque	Angle	Offset	Tube	Reference Code		Reference Code	
					Slot 0,18		Slot 0,22	
					Right	Left	Right	Left
Edgewise Standard	0	0	0	3.6	353-007	353-008	353-107	353-108
MBT 1° Bicuspid	-20°	0	0	3.6	353-295	353-296	353-395	353-396
MBT 2° Bicuspid	-10°	0	0	3.6	353-293	353-294	353-393	353-394
Roth	-25°	5	0	3.6	353-233	353-234	353-335	353-336

## 12.4. Weldable Buccal Tube

### UPPER Tooth

	Prescription	Torque	Angle	Distal Offset	Length (mm)	Occlusal or Gingival	Hdgr Tube	Reference Code			
								Slot 0,18		Slot 0,22	
								Right	Left	Right	Left
 1° Molar Upper- Triple Convertible with Hook Ricketts NO Hook Univ.	Roth	-10	0	10	3.5	O	.045	340-830	340-831	340-919	340-920
	MBT	-14	0	10	3.5	O	.045	340-865	340-866	340-965	340-966
	Ricketts	0	0	0	3.5	G	.045	340-011	340-012		
	Ricketts	0	0	0	3.5	G	.045	340-013	340-013	-	-
 1° Molar Upper - Double Convertible whit Hook	Roth	-10	0	10	3.5	O	.045	326-825	326-826	326-925	326-926
	Roth	-10	0	10	3.5	G	.045	326-829	326-830	326-931	326-932
	MBT	-14	0	10	3.5	O	.045	326-871	326-872	326-971	326-972
	MBT	-14	0	10	3.5	G	.045	326-859	326-860	326-892	326-893
	Edge STD	0	0	0	3.5	O	.045	326-805	326-806	326-905	326-906
 Ricketts Not Convertible	Ricketts	0	0	0	3.6	O	.045	342-011	342-011	-	-
 1° Molar Upper- Single Convertible Whit Hook	Roth	-10	0	10	3.5	-	-	336-877	336-878	336-977	336-978
	MBT	-14	0	10	3.5	-	-	336-873	336-874	336-973	336-974
	Edge STD	0	0	0	4.0	-	-	300-201	300-201	300-301	300-301
 2° Molar Upper- Single Not convertible with Hook	Roth	-10	0	10	3.5	-	-	336-877	336-878	336-977	336-978
	MBT	-14	0	10	3.5	-	-	336-873	336-874	336-973	336-974
	Edge STD	0	0	0	3.5	-	-	300-201	300-201	300-301	300-301

### LOWER Tooth

	Prescription	Torque	Angle	Distal Offset	Length (mm)	Occlusal or Gingival	Hdgr Tube	Reference Code			
								Slot 0,18		Slot 0,22	
								Right	Left	Right	Left
 1° Molar Lower Double Rect.	Roth-Rect.	-25	0	5°	3.5	G	-	342-813	342-814	342-913	342-914
	MBT-Rect.	-20	0	0	3.5	G	-	342-815	342-816	342-915	342-916
 1° Molar Lower- Double Round Conv. whit Hook	Roth-Round	-25	0	5°	3.5	G	.045	339-813	339-814	339-913	339-914
	MBT-Round	-20	0	0	3.5	G	.045	339-815	339-816	339-915	339-916
	Edge STD	0	0	0	3.5	G	.045	339-803	339-804	339-903	339-904
 Ricketts Not Convertible	Ricketts	0	0	0	3.6	O	.045	342-011	342-011	-	-
 1° Molar Upper- single Convertible Whit Hook	Roth	-25	0	5	3.5	-	-	336-813	336-814	336-913	336-914
	MBT	-20	0	0	3.5	-	-	336-815	336-816	336-915	336-916
 2° Molar Upper- Single Not convertible with Hook	Roth	-25	0	5°	3.5	-	-	324-231	324-232	324-335	324-336
	MBT	-10	0	0	3.5	-	-	324-293	324-294	324-393	324-394
	Edge STD	0	0	0	3.5	-	-	324-001	324-002	324-101	324-102