

Sterilization Instructions

All Instruments must be cleaned and sterilized in the open position.

Never exceed the temperature of 246°C (475°F).

Sterilize for 5 minutes at 134°C (270°F) or for 20 minutes at of 121 °C (250°F)

Lubricate Joints.

Warranty

We Warranty every instrument with our brand name to be free of functional defects in workmanship and materials when used normally for its intended use. Any of our products that prove to be defective in material and workmanship will be repaired or replaced at our discretion. This guarantee (warranty) does not cover damage caused by abuse or misuse. Guarantee is voided if unauthorized service resources repair or modify the instruments.