

30.0. Clear Retainer after Braces



- Superior Aesthetics
- Unparalleled Durability
- Gentle on teeth
- BPA Free
- Virtually unbreakable
- Virtually resistant to stains
- Unaffected by sunlight



The ClearBow is a virtually transparent, encouraging patient to wear their orthodontic retainers and show off their beautiful new smile without feeling self-conscious.

Now you can confidently show off your beautiful smile with the durable ClearBow



Hawley retainer; a new clear retainer that is fabricated with a virtually unbreakable, BPA-free plastic orthodontic wire (ClearBow). The ClearBow orthodontic wire is made in Canada and is not damaged by boiling water or stained by coffee, cola, beet juice or red wine.



The ClearBow Hawley retainer is designed to provide superior aesthetics encouraging patients to wear their durable clear Hawley retainer without feeling self-conscious.

Unlike other types of clear retainers, the ClearBow Hawley retainer covers little or no biting surfaces of the teeth. When properly designed by the orthodontist, the

ClearBow Hawley retainer allows the teeth to settle into occlusion, which will also permit even distribution of the biting forces.

The ClearBow orthodontic wire is guaranteed not to break when the ClearBow Hawley retainer is used under normal conditions specified by the doctor. The ClearBow Hawley retainer can be prescribed to children, adolescents and adults by a general dentist or orthodontist.

AE500-201	ClearBow non-extraction (Size +2) Pkg.5	AE500-205	ClearBow extraction (Size -2) Pkg.3
AE500-202	ClearBow non-extraction (Size +3) Pkg.5	AE500-206	ClearBow extraction (Size -3) Pkg.3
AE500-203	ClearBow non-extraction (Size +4) Pkg.5	AE500-207	ClearBow extraction (Size -4) Pkg.3
AE500-204	ClearBow non-extraction (Size +5) Pkg.5	AE500-208	ClearBow extraction (Size -5) Pkg.3
AE500-210	ClearBow Kit Extraction Pkg 8	AE500-211	ClearBow Intro Kit Pkg 24
AE500-209	ClearBow Kit Non-Extraction Pkg 8	AE500-212	ClearBow Intro Kit Instrument

